


S4PE TEST SCRAPERS

Part # M266 EA5 000

The test scraper, mounted on the S4PE, enables an operational check (bump test) of the AP2C and AP4C.



It is impregnated with a substance which, when heated and vaporized by the S4PE activates AP2C and AP4C's display indicators (light up on both Phosphorus and Sulfur channels).

• Packaging

The scrapers are supplied in a box of 10. They are protected with a cover. They can be opened and the scrapers used one by one. The boxes of scrapers can be stored unopened for 10 years. A scraper can be used up to 10 times.

They differ from the sampling scrapers because they are blue (sampling scrapers are black).

• Characteristics

- Size of a test scraper (LxWxH): 12,3 mm x 10 mm x 18 mm (0.48" x 0.39" x 0.7")
- Size of a box of scrapers (LxWxH): 150 mm x 17 mm x 26 mm (5.9 " x 0.66" x 1.02")
- Weight of one scraper: 1 gr (0.03 oz)
- Weight of full box of scrapers: 32 gr (1.1 oz)
- Operating temperature: -32°C to +55°C (-25.6°F to +131°F)
- Storage temperature: -39°C to +71°C (-38.2°F to +160°F)

